
Model	A710-2
Rated voltage	110-240V~ 50/60Hz
Rated power	220W
Rated time	20min
G.W	132.2KG
Noise value	≪60dB
Safety structure	Class I






# A710-2 MASSAGE CHAIR



# 23 types of automatic massage programs



### Al massage robots



135cm
"S+L"shape long rail

**1.5-21**cm width

12cm

n dej





### Shawl hot compress

It is provided with a unique shawl assembly, with innovative heating functions on the shoulder, back, waist and abdomen and heating with graphene (a wormwood herb medicine bag is placed at the bottom of the shawl). Users can place the shawl at different massage positions according to massage requirements.







wormwood herb medicine bag



asseml

# Leg and foot massage

#### calf kneading function

kneading and airbag squeezing around the inside and outside of the calves

### foot roller massage

scrape and massage mastoid at the arch and toes

#### ankle stretching function

use the airbag squeezing function, stretch the ankle by stretching the foot rest up and down

## heating function at the soles of the feet and toes

use the airbag squeezing function, stretch the ankle by stretching the foot rest up and down

